



IDENTIFICATION OF GY SERIES

■ INSERT




■ CBN INSERT




★1 Select a seat size with the same symbol as that of modular blade and mono block holder.

■ MODULAR BLADE

● EXTERNAL and INTERNAL



● FACE GROOVING



★1 Select a seat size with the same symbol as that of the insert.

★2 The maximum groove depth is a value when used for external grooving and changes according to the insert used.
For internal grooving, refer to the maximum groove depth on pages F112 – F116.

★3 GYM20R/LA-○10, GYM20R/LA-○12, GYM25R/LA-○12 and GYM25R/LA-○14 can be used for both external and internal grooving.